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INDICE / SOMMAIRE

INSTALATION

Requirements for optimum installation

REQUIREMENTS FOR AN OPTIMAL INSTALLATION

To obtain the maximal performance and the best results, we recommend installing the balance according to the following requirements:

- a) Install it on a firm and stable surface.
- b) Place it on a stable base free of vibrations.
- c) Don't expose it to direct sunlight.
- d) Don't submit it to atmospheres with corrosive gases.
- e) Keep it in an environment free of dust and dirt.
- f) Work at temperatures between 15° and 30° C.
- g) Use the balance in environments with humidity lower than 95% HR.
- h) Keep the balance away from equipment that produces magnetic and acoustical interferences
- i) Calibrate it in its definitive location and check it regularly with the optional calibration weight, if necessary, recalibrate it.

Unpacking

Unpack the balance carefully. Inside the carton you will find the following components:

1. This operation manual (please, read it before to use the balance).
2. A balance, series STZ
3. A stainless steel pan
4. An AC/Adaptor of 12 V / 800 mA; don't use any other adaptor, it could damage the balance.

VERY IMPORTANT: Keep the original package. It could be useful in the future in transporting the balance.

INSTALLATION

Operating temperature: 5°C – 35°C

Range of temperature fluctuation: $\leq 2^{\circ}\text{C}/\text{h}$

Relative humidity: 50%-85%

1. Place the balance on a firm and stable surface.
2. Connect the AC/Adaptor
3. Press the  key to turn on the balance, the display will show "8.8.8.8.8.8", the maximal capacity, "S—CAL" and finally "0.0" or "0.00"

4. To obtain maximal accuracy, before using the balance, we recommend keeping it turned on for 10 minutes.

KEYBOARD DESCRIPTION

 To turn on or turn off the balance.

 Tara key

 Functions key

 Print key

 Backlight key

 Calibration key

WEIGHT UNIT SELECTION

1. Keep the **Mode** key pressed for 3 minutes.
2. The message “**Unit**” will light.
3. Press the **→T←** key to confirm.
4. The message “**Unit**” will light again.
5. Press the **Mode** key to choose the desired weight unit.
6. To confirm, press the **→T←** key.

TARE

This function allows obtaining the net weight of a product.

1. With display showing zero value, place a container on the weighing pan.
2. Press the **→T←** key.
3. The balance will subtract the container’s weight, put the product into the container and then read the product’s weight on the display (net weight).

To return to normal weighing mode, remove the container from platter and press the **→T←** key.

PIECES COUNTING FUNCTION

1. Press the **Mode** key for 3 seconds.
2. The display will show the message “Unit”
3. Press again the **Mode** key.
4. The display will show the message “Count”
5. Press the **→T←** to confirm.
6. The message “pcs” will appear on the display.
7. Press the **Mode** key.
8. Select the sample number from “10pcs” to “500pcs”.
9. Put on the platter the quantity of pieces that corresponds to a selected value.
10. Press the **→T←** key.
11. The display will show “- - - - -” for a few seconds.
12. The balance is ready to perform the counting of pieces.

The greater quantity of pieces put in the initial sample, the greater accuracy in balance's results.

If the sample value is less than 0.2d, the display will show the message **"no-Cou"**. This means that the sample is too small; therefore, you will need to add a greater sample.

PERCENTAGES FUNCTION (%)

1. Keep the  key pressed for 3 seconds.
2. Press again the  key twice.
3. The message "PER" will appear on the display.
4. Press the  key to confirm.
5. The display will show "100%"
6. Put on the platter an object considered as 100%.
7. Press the  key
8. Remove the object from the weighing pan.
9. Place other items and the balance will calculate the percentage in regard to the default sample.

If the sample value is less than 0.2d, the display will show the message

"no-PER". This means that the sample is too small; therefore, you will need to add a greater sample.

BACKLIGHT

Press the  key to activate or deactivate the internal backlight.

LOW BATTERY INDICATOR

If the battery indicator is activated in the upper left side on the display, this means that the battery is low. In this case, recharge the battery using the provided adaptor.

SPEED SELECTION

Press the  key for 3 seconds to select the desired mode: SPEED 1 or SPEED 2.

AUTO-ZERO FUNCTION

1. When the balance is Off, keep pressed the  key and without releasing it, press the  key until the screen shows the message -Zero-, now release both keys.
2. Press and hold the  key until the display shows the message “Zero*d”
3. Press again the  key.
4. Select the desired option (from 0 to 5)
 - Zero0d Auto-zero function not activated.
 - Zero5d Auto-zero function activated = 0.05g – 300.09g
0.05g – 600.09g
0.05g – 3000.09g

(Default value: Zero5d)

AUTO-TARE FUNCTION

1. Press the  key, the display will show the message “-tArE-“
2. Press the  key, the display will show the message “tArE*d”
3. Press the  key twice.
4. Select the desired option (from 0 to 9)
 - a. tArE0d - Auto-tare function not activated.
 - b. tArE9d - Auto-tare function activated. This means that when a weight is placed on the weighing pan, the balance will tare automatically its weight.

(Default value: tArE9d)

5. Press the  key to confirm and save the configuration. The balance will return to normal weighing mode “0.00”

CALIBRATION

1. Press the  key for 3 seconds.
2. The display will show “d3”
3. If your balance is the AHZ-300 model, press the  key until the screen shows the message “d2”, then, the display will show “200g”. Put on the weighing pan a calibration weight of 200g.
4. If your balance is the AHZ-600 model, press the  key until the screen shows the message “d5”, then, the display will show “500g”. Put on the weighing pan a calibration weight of 500 g.
5. If your balance is the AHZ-3000 model, press the  key until the screen shows the message “d3”, the display will show “3000g”. Put on the weighing pan a calibration weight of 3000 g.
6. The display will show automatically “-CAL-“
7. The display’s reading will remain with the value of calibration weight.
8. The calibration procedure has finished.

LINEARITY CALIBRATION

1. Press the  key for 3 seconds.
2. Press again the  key until the screen shows the message “CAL---L”
3. If your balance is the AHZ-300 model, the display will show “100g”. Put on the weighing pan a calibration weight of 100g.
4. If your balance is the AHZ-600 model, the display will show “200g”. Put on the weighing pan a calibration weight of 200 g.
5. If your balance is the AHZ-3000 model, the display will show “1000g”. Put on the weighing pan a calibration weight of 1000 g.
6. The display will show automatically “-CAL-“
7. The display’s reading will remain with the value of calibration weight.
8. The calibration procedure has finished.

PRINTER FUNCTION

By using the RS-232C output, connect the optional printer.

Press the  key to print displayed data.

PRINT CONFIGURATION

Mode 1 “hAnd” - By pressing the  key, the printer prints the transmitted data

Mode 2 “Auto” – The data will be transmitted automatically when the weight placed on the pan is > 3d and the display reading is stable.

Mode 3 “Contin” – The data is continuously transmitted to the printer.

1. Press the  key for 3 seconds.
2. The display will show the message “Unit”
3. Press the  three times.
4. The display will show the message “Prt”
5. Press the  to confirm.
6. The display will show the message “hAnd”
7. Press the  key to select the desired mode: “hAnd”, “Auto” and “Contin”
8. Press the  key to confirm the selected mode.

* Parameter of RS-232 interface:

Format of data:

Each 10BIT as a data, the format:

BIT1, BIT2, BIT3,,BIT8, BIT9, BIT10

BIT1: start bit

BIT2----BIT9: data bit

BIT10: final bit

Baud rate: 9600 bps

GUARANTEE

This balance is guaranteed for one year from the delivery date. The guarantee covers any fabrication defect of the material.

During this period, GRAM PRECISION covers the manpower and spare parts costs for the reparation of the scale.

This guarantee does not cover the failures caused by an inappropriate use or overload.

The guarantee does not cover the freight cost (transport) necessary to repair the balance.

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